



Nihon Kohden America, Inc. announces the auction of its latest EEG technology to support the AES Lennox and Lombroso Trust Fund/Susan Spencer Research Fund
EEG-1200A Diagnostic and Monitoring Solutions



System includes cart, camera and microphone.

100% of proceeds will be presented to the AES at the annual meeting in San Antonio to support the AES Lennox and Lombroso Trust Fund/Susan Spencer Research Fund.

Link to Lennox and Lombroso Trust Fund/Susan Spencer Research Fund:
www.aesnet.org/index.cfm?objectid=A8B423A2-E7FF-0F41-28973682F5AC6369

A \$38,250 value, the EEG-1200 auction is in effect from September 1 through November 1, 2010.

Minimal bid is \$10,000.00

Download the bid document from www.nkusa.com/neurology_cardiology/
Insert your name and contact information along with your bid amount and email to AES_EEGAuction@nkusa.com.

Delivery and bid payment to be completed by November 30, 2010.

The highest bid amount will be posted every Monday on the Nihon Kohden America web site (above) through October 25th.

Please contact Kathy Hart, Director of Neurology at NKA with questions or more information at 1-800-325-0283, ext 3405 or email at Kathy_hart@nkusa.com

See reverse page for the EEG-1200 auction configuration and options. If options are desired, contact NKA for a separate pricing proposal. Normal purchase terms and conditions apply.

Auction and purchase offered only in the United States.

The following information is required:

QUOTATION # KH100616-01

Name	Company Name	
Address		Terms: 100% Due upon Delivery
		Warranty: 1 Year Parts & Labor
City, State & Zip Code		FOB: Destination
Phone #	Fax #	Delivery by November 30, 2010
Email		AES_EEGAuction@nkusa.com

QTY	PART NUMBER	AES AUCTION DESCRIPTION - \$38,250 VALUE Acquisition Systems
1	EEG-1200A	Neurofax EEG-1200A Acquisition system. Includes PC (SFF) - Intel Core 2 Duo, 3.0 GHz, 1333 FSB, 2GB Ram, 250 GB Hard Drive, DVD +/- RW, CD-R, Keyboard & Mouse, Gigabit Ethernet NIC, Patient Isolation, Microsoft Windows XP Professional with NeuroWorkbench, EEG-1200 Acquisition, Review, File Utility software.
1	Choice of Carts	NK-NIH8000 or KE-122A System Cart
1	NK-FLTSCREEN-19	19" Flat Screen TFT Monitor (1280 x 1024)
1	NK-AUDIOPACK	Amplifier and Speakers
1	NK-251A	Sony IP Camera Mount for KE-122A and Modo carts
1	JE-921A	Amplifier 32 Channel w/ SP02 and CO2 Capability
1	JE-922A	Mini input box for JE-921A, includes 5m cable
1	NK-POUCH-JE922	Mini Input Box Pouch for JE-116, JE-113 and JE-914, includes shoulder strap
1	NK-MPEG4-SW	Digital Video Software with MPEG4 32, includes QI-120A USB Video Capture
1	NK-CAMERA-1	SONY IP camera (Pan, Zoom, Tilt)
1	NK-MICROPHONE2	Microphone system
1	LS-120AJ	EEG-1200 Photic Stimulation Control Unit
1	LS-703A	Photic Flash Lamp Assembly, EEG-1200, EEG-9000 only
1	NK-CDNOW	(Port-a-View) Records EEG and Video to CD for viewing on Windows PC
1	NK-Installation	Installation and Termination per unit up to 20 units
1	NK-EEG-DV-Training	Up to 4 days including travel days and configuration
1	NK-BMET	Tuition only, Biomedical System Training

BID AMOUNT \$

This Quotation is Valid for a period of 60 days.
Available only in the United States.

OPTIONS NOT INCLUDED IN PRICE - PLEASE CONTACT NKA FOR OPTION PRICING

QTY	PART NUMBER	DESCRIPTION	UNIT PRICE	TOTAL PRICE
1	EEG-1200R	Neurofax EEG Review system. Includes PC (SFF) - Intel Core 2 Duo, 3.0 GHz, 1333 FSB, 2GB Ram, 250 GB Hard Drive, DVD +/- RW, CD-R, Keyboard & Mouse, Gigabit Ethernet NIC, Microsoft Windows XP Professional with Review, NeuroWorkbench, NWB Report Generation and File Utility software.		
1	NK-NEUROFAX-R	Neurofax EEG Review system. Includes PC (Tower - Dual Monitor Support) - Intel Core 2 Duo, 3.0 GHz, 1333 FSB, 2GB Ram, 250 GB Hard Drive, DVD +/- RW, CD-R, Keyboard & Mouse, Gigabit Ethernet NIC, Microsoft Windows XP Professional with Review, NeuroWorkbench, NWB Report Generation and File Utility software.		
1	NK-FS-19-HR	20" Flat Screen Hi Resolution TFT Monitor (1600 x 1200)		
1	NK-FLTSCREEN-24	24" Wide Flat Screen Hi-Res. TFT Monitor (1920 x 1200)		
1	NK-REVIEW-SITE	EEG-1100 Review Software Site License (Single User)		
1	NK-CPR	Persyst Comprehensive EEG Review Station Pack		
1	NK-CPA	Persyst Comprehensive EEG Monitoring Pack		



90 Icon Street, Foothill Ranch, CA 92610

(800) 325-0283 • (949) 580-1555 • Fax: (949) 580-1550 • www.nkusa.com