

# SECTION 1

## Electrode Application Supplies



Adhesive

Skin Preparation Cream

Gel

Accessory Order Form

## Electrode Application Supplies

### Elefix, Hard Tub

Electrode adhesive, 400 grams per tub. Used with all surface applied electrodes. Shelf life is 3 years unopened; 6 months after opening.

Applicable Instruments: All EEG and EP/EMG Instruments

**Part No: Z-401CE**

**Qty: 3/BOX**



### Elefix, Squeeze Tube

Electrode adhesive, 180 grams per tube. Used with all surface applied electrodes. Shelf life is 3 years unopened; 6 months after opening.

Applicable Instruments: All EEG and EP/EMG Instruments

**Part No: Z-181JE**

**Qty: 10/BOX**



### SkinPure

Skin preparation cream, 135 grams per tube. Used with all surface applied electrodes. Shelf life is 3 years unopened; 6 months after opening.

Applicable Instruments: All EEG and EP/EMG Instruments

**Part No: YZ-0019**

**Qty: 10/BOX**



### Aquasonic 100

Ultrasonic transmission gel, 8.5 ounce tube. Use for all EEG or EMG procedures where a viscous gel is required. Completely aqueous, will not stain clothing.

**Part No: KK-01-08**

**Qty: 1**

### SignaGel

Multipurpose saline electrode EMG gel, 8.0 ounce tube

**Part No: KK-15-25**

**Qty: 1**



## Section

1

Electrode Application Supplies

2

EEG/PSG

3

EP/EMG/NCV

