

Contact:

Jeannie Callahan
Nihon Kohden America, Inc.
jeannie_callahan@nkusa.com
(949) 580-1555 ext.3374
www.nkusa.com

FOR IMMEDIATE RELEASE**Introducing the QP-160 cEEG Trending Software**

Foothill Ranch, CA December 10, 2010 – Nihon Kohden today announced the immediate release of the QP-160 cEEG trending software for the EEG-1200 system.

The use of continuous EEG (cEEG) can enhance the neurological assessment and care of critically ill patients through early detection of reversible neurological deterioration. The Nihon Kohden superior signal quality amplifiers ensure reliable cEEG trends. The QP-160 software is designed to fully integrate the EEG data with user selectable trends over time. Areas of clinical interest are easily recognized at the bedside and provide clinicians fast review capabilities of the associated EEG data over days, and even weeks, in a single display with synchronized audio and video data.

High Quality Trends included in this application are DSA (Density Spectral Array), FFT (Fast Fourier Transform), BSR (Burst Suppression Ratio), Alpha/Delta Plot (FFT Power Ratio), aEEG (Amplitude Integrated EEG), IBI (Inter Burst Interval) and BPM (Burst per Minute). In addition the QP-160 provides important values for DC signals such as SaO2 and etCO2. Trends can be created in both a Referential and Bipolar configuration and all are integrated within both Nihon Kohden Acquisition and Review applications.

The versatile, comprehensive QP-160 application is easy to use and offers:

- Integration with the EEG-1200 system acquisition and review.
- Nihon Kohden's superior signal quality amplifiers to ensure reliable cEEG trends.
- Ability to view real time EEG and a variety of single or dual time base trend displays remotely.
- Remote viewing and reviewing of the entire long term recording.
- Benefits to the neurologist and neuro intensivist for rapid assessment and intervention.

About Nihon Kohden America

Nihon Kohden America was founded in 1979 as a wholly owned subsidiary of Nihon Kohden Corporation, Japan (founded in 1951). Nihon Kohden is a market leader in Patient Monitoring, Sleep Assessment, Neurology and Cardiology instrumentation. Our products are distributed to hospitals, physician's practices, mobile services, and outpatient and surgery centers throughout North America. More information about Nihon Kohden can be found at www.nkusa.com.

#####