

SECTION 3

EP/EMG/NCV



Electrodes

Needle Electrodes

Injectable Needle Electrodes

Stimulators

Cables/Connectors

Miscellaneous


Accessory Order Form

Reusable EMG/NCV Electrodes

Adult Finger Clip Electrode - 24" red/black lead, stainless steel clips

Part No: KK-FC-EZ

Qty: 1



Single Alligator Clip (Black) - 24" lead

Part No: KK-ETL-24BSAF

Qty: 1



Single Alligator Clip (Green) – 24" lead

Part No: KK-ETL-24GSAF

Qty: 1

Single Alligator Clip (Red) - 24" lead

Part No: KK-ETL-24RSAF

Qty: 1

Braided Alligator Clips – R, G, BK; 100cm lead

Part No: KK-ETL-RGB1

Qty: 1

Ring Electrode Pair – Neoprene slide adjustments, stainless steel springs, 24" red/black lead

Part No: KK-RE(R)-SAF


Qty: 1



Ground Electrode – Nickel-plated brass, 32mm Flat Disc, 24" lead

Part No: KK-285-32-24TP

Qty: 1



Ground Electrode – Nickel-plated brass, 32mm Flat Disc, 48" lead

Part No: KK-285-32-48TP

Qty: 1

Section

①

Electrode
Application
Supplies

②

EEG/PSG

③

EP/EMG/NCV

Paired Disc Electrode – 9mm round, locator pin, surgical steel, 24” red/black lead

Part No: KK-SEF-924SAF

Qty: 1



SlimLine Bar Electrode – Flat 9mm discs, 30mm spacing, 24” red/black lead

Part No: KK-DDB-F30

Qty: 1



Reusable EEG Electrode, Gold, 9mm

Part No: KK-WBT-9GSAF

Qty: 10/PKG



EMG/NCV Monopolar Start-Up Kit

Includes: KK-DDBF-30, KK-RE(R)-SAF, KK-SEF-924SAF, KK-ETL-RGB1, KK-ETL-24BSAF, KK-ETL-24GSAF, KK-ETL-24RSAF, KK-285-32-24TP, KK-PT-30, NK-RE-D, KK-CA-510, KK-CA-610, KK-715-08-K/C (X3), Measuring Tape, Variety of Monopolar Samples, and Aquasonic 100, 8.5oz Tube

Part No: NK-STARTUP-MONO

Qty: 1

EMG/NCV Concentric Start-Up Kit

Includes: KK-DDBF-30, KK-RE(R)-SAF, KK-SEF-924SAF, KK-ETL-RGB1, KK-ETL-24BSAF, KK-ETL-24GSAF, KK-ETL-24RSAF, KK-285-32-24TP, KK-PT-30, KK-RE-D, KK-CA-510, KK-CA-610, KK-715-08-K/C (X3), KK-1741, Measuring Tape, Variety of Concentric Samples, and Aquasonic 100, 8.5oz Tube

Part No: NK-STARTUP-CONC

Qty: 1

Electrode Set, Shielded EP

Shielded electrodes with 60-inch lead cable and Touch Proof connector; Package includes one NE-432B jumper cable

Applicable Instruments: MEB-4200A, MEB-2200A, MEB-5500, MEB- 9100, MEB-9200/9300, MEB-9400

Part No: NE-132B

Qty: 5/PKG



Section

1

Electrode
Application
Supplies

2

EEG/PSG

3

EP/EMG/NCV

Disposable EMG/NCV Electrodes (one time use only)

Tab Electrodes, 15mm x 20mm, 150cm lead; requires reusable alligator clips

Part No: KK-700-15-K/C

Qty: 12/PKG



| Part Number | Description | Qty |
|---------------|--|---------|
| KK-PT-30 | Bar Electrodes - 10mm x 10mm, requires reusable alligator clips | 100/BOX |
| KK-714-15/M1 | Large Ground Electrode | 1 |
| KK-RE-D | Ring Electrodes - 8mm x 102mm, requires reusable alligator clips | 100/BOX |
| KK-715-05-K/C | Small Sensor NCS, 50cm disposable lead wire | 12/BOX |
| KK-CA-310 | Ground Electrodes - 28mm x 20mm, requires reusable alligator clips | 100/BOX |
| KK-715-08-K/C | Small Sensor NCS, 80cm disposable lead wire | 12/BOX |
| KK-CA-510 | Ground Electrodes, 34mm x 28mm, requires reusable alligator clips | 100/BOX |
| KK-715-12-K/C | Small Sensor NCS, 120cm disposable lead wire | 12/BOX |
| KK-CA-610 | NCV Electrodes, Round, 10mm, requires reusable alligator clips | 100/BOX |


Subdermal Needle Electrodes (single use only)

| Part Number | Description | Qty |
|-------------|--------------------------|-------------|
| KK-S43-438 | 13mm x 27ga (100cm lead) | Qty: 25/PKG |
| KK-S43-638 | 13mm x 27ga (150cm lead) | Qty: 25/PKG |
| KK-S43-838 | 13mm x 27ga (200cm lead) | Qty: 25/PKG |
| KK-S43-938 | 13mm x 27ga (250cm lead) | Qty: 25/PKG |

Twisted Pair 13mm x 27ga, 150cm lead

Part No: KK-S44-637

Qty: 40/PKG



Twisted Pair 13mm x 27ga, 200cm lead

Part No: KK-S44-837

Qty: 40/PKG

Corkscrew, 100cm lead

Part No: KK-001-220

Qty: 6/PKG



Section

1

Electrode
Application
Supplies


2

EEG/PSG

3

EP/EMG/NCV

Concentric EMG Needles




| Part Number | Description | Qty |
|--------------|-------------|--------|
| KK-740-25-30 | 25mm x 30ga | 25/PKG |
| KK-740-38-45 | 38mm x 26ga | 25/PKG |
| KK-740-25-45 | 25mm x 26ga | 25/PKG |
| KK-740-50-45 | 50mm x 26ga | 25/PKG |
| KK-740-30-35 | 30mm x 28ga | 25/PKG |
| KK-740-75-65 | 75mm x 23ga | 25/PKG |

Concentric Needle Cable

Reusable Concentric EMG Needle Cable, 5 Pin DIN, 1m

Part No: **KK-1741**

Qty: 1



Monopolar EMG Needles

| Part Number | Description | Qty |
|--------------|------------------------------|--------|
| KK-742-25-36 | Detachable Lead, 25mm x 28ga | 40/PKG |
| KK-743-25-36 | Lead Attached, 25mm x 28ga | 40/PKG |
| KK-742-38-36 | Detachable Lead, 38mm x 28ga | 40/PKG |
| KK-743-38-36 | Lead Attached, 38mm x 28ga | 40/PKG |
| KK-742-38-45 | Detachable Lead, 38mm x 26ga | 40/PKG |
| KK-743-38-45 | Lead Attached, 38mm x 26ga | 40/PKG |
| KK-742-50-45 | Detachable Lead, 50mm x 26ga | 40/PKG |
| KK-743-50-45 | Lead Attached, 50mm x 26ga | 40/PKG |
| KK-742-75-50 | Detachable Lead, 75mm x 25ga | 40/PKG |
| KK-743-75-50 | Lead Attached, 75mm x 25ga | 40/PKG |

Monopolar Needle Cable

Reusable Monopolar Needle Cable, 61cm length

Part No: **KK-1744**

Qty: 1

Section

1

Electrode
Application
Supplies


2

EEG/PSG

3

EP/EMG/NCV

Injectable Monopolar EMG Needles



Ambu Neuroline® Inoject™ Injectable Monopolar EMG Needles (single use only)

| Part Number | Description | Qty |
|--------------|------------------------|--------|
| KK-744-25-30 | Inoject™ - 25mm x 30ga | 10/PKG |
| KK-744-30-36 | Inoject™ - 30mm x 28ga | 10/PKG |
| KK-744-35-40 | Inoject™ - 35mm x 27ga | 10/PKG |
| KK-744-38-45 | Inoject™ - 38mm x 26ga | 10/PKG |
| KK-744-50-50 | Inoject™ - 50mm x 25ga | 10/PKG |
| KK-744-75-55 | Inoject™ - 75mm x 24ga | 10/PKG |

Chalgren Injectable Monopolar EMG Needles (single use only)



| Part Number | Description | Qty |
|--------------|------------------------|--------|
| KK-242-025TP | Chalgren - 25mm x 30ga | 10/PKG |
| KK-242-375TP | Chalgren - 75mm x 23ga | 10/PKG |
| KK-242-550TP | Chalgren - 50mm x 25ga | 10/PKG |
| KK-242-637TP | Chalgren - 37mm x 26ga | 10/PKG |
| KK-242-725TP | Chalgren - 25mm x 27ga | 10/PKG |

Section

1

Electrode
Application
Supplies

2

EEG/PSG

3

EP/EMG/NCV

Stimulators

Disposable Handle for Pedicle Hole Stimulator

Part No: KK-RMA1012


Qty: 5/PKG



Reusable Pedicle Hole Stimulator 120cm Lead

Part No: KK-RMA1009

Qty: 1



Disposable Pedicle Hole Stimulator

Part No: KK- RMA1009D

Qty: 12/PKG

EMG Cables

| Part Number | Description | Qty |
|-------------|----------------------------------|-------|
| 584351 | USB Cable for MEB-9100 System | 1 |
| AAA-15482 | Cable for Input Box JB-281BK, 9M | 1 |
| AAA-15955 | USB Cable for MEB-92/93/2200, 9M | 1 |
| BM-001B | 5 Pin MEB Cable | 1 |
| C927 | Leg Wire Extensions; 48-inch | 5/PKG |

Cable, Stimulation/Input

5-foot extension cable, 5 Pin, round DIN connector allows for stimulator connection or electrode input connection; 1.5mm DIN type Touch Proof lead input block provides for simultaneous negative and positive input connection

Applicable instruments: MEB and MEM Series Instruments

Part No: SP-166-5

Qty: 1

25-foot extension cable. 5 Pin, round DIN connector allows for stimulator connection or electrode input connection; 1.5mm DIN type Touch Proof lead input block provides for simultaneous negative and positive input connection

Applicable instruments: MEB and MEM Series Instruments

Part No: SP-166-25

Qty: 1

5-foot extension cable. 8 Pin, round DIN connector allows for stimulator connection or electrode input connection; 1.5mm DIN type Touch Proof lead input block provides for simultaneous negative and positive input connection

Applicable instruments: MEB and MEM Series Instruments

Part No: SP-168-5

Qty: 1

Section

1

Electrode
Application
Supplies

2

EEG/PSG

3

EP/EMG/NCV

25-foot extension cable; 8 Pin, round DIN connector allows for stimulator connection or electrode input connection; 1.5mm DIN type Touch Proof lead input block provides for simultaneous negative and positive input connection

Applicable instruments: MEB and MEM Series Instruments

Part No: SP-168-25

Qty: 1

Cable, Adaptor, Extension

11-inch extension cable for connection to the 5 pin DIN connections on all Neuropack input boxes; allows for simultaneous positive, negative, and ground connections which accept Touch Proof leads

Applicable instruments: MEB and MEM Series Instruments

Part No: NK-TPCAB5PIN

Qty: 1

Cable, BNC to BNC

3-foot cable with BNC to BNC male connectors; connects VM-125 or PM-155 Visual Monitors to all Nihon Kohden Neuropack products

Applicable instruments: MEB and MEM Series Instruments

Part No: 5050-B-60-0

Qty: 1

Cable, Gelpad

36-inch Twisted-Pair Alligator Clip leads with 1.55mm DIN type connectors. Use with NK-GELPAD pre-gelled electrode pads

Applicable instruments: MEB and MEM Series Instruments

Part No: NK-GELPAD-TPCBL

Qty: 1

Miscellaneous

Buzz Kill Noise Detector

OR Version

Part No: KK-BK-OR

Qty: 1



Disposable Temperature Strips

Part No: KK-DTS01

Qty: 50/PKG

Cable, Jumper

DIN Shortening lead cable for NE-132B Electrode Set

Part No: NE-432B

Qty: 5/PKG

Section

1

Electrode
Application
Supplies

2

EEG/PSG

3

EP/EMG/NCV


Felt Tip

Replacement felt tips for use with NM-420S bar stimulator

Applicable instruments: All MEB and MEM Series Instruments

Part No: 5010097

Qty: 100/PKG



Holder, Input Box

Quick disconnect bracket for JB-482BK, JB-216BK, JB-241BK, JB-281BK, MS-210BK input boxes. Includes one 6143-001909 holder assembly.

Applicable Instruments: MEB-4200, MEB-5508A, MEB-2200, MEB-9200/9300

Part No: DI-0016

Qty: 1

Ear Phone

4 Pin Connector Earphone

Applicable Instruments: MEB-4200A, MEB-5304K, MEB-7102A, MEM-7102A, MEM-7102A, MEM-4104K, MEB-2200

Part No: YE-102J

Qty: 1

Ear Phone

Applicable Instruments: MEB-2200A, MEB-9200/9300, MEB-9400

Part No: YE-103J

Qty: 1

Foot Switch

Application Instruments: MEB-9100 MEB-9102, MEB-9104, MEB-9200/9300

Part No: RY-480B

Qty: 1

Goggles, LED

Applicable Instruments: MEB-4200A, MEB-5304K, MEB-7102A, MEM-4104K, SLS-3500, MEB-2200

Part No: LS-101J

Qty: 1

Goggles, LED

Application Instruments: MEB-9100, MEB-9200/9300, MEB-9400

Part No: LS-102J

Qty: 1

Section

1

Electrode
Application
Supplies

2

EEG/PSG

3

EP/EMG/NCV

Headphones

4 Pin Connector Headphones

Applicable Instruments: MEB-4200A, MEB-5304K, MEB-7102A, MEM-7102A, MEM-4104K, MEB-2200

Part No: DR-531B-10

Qty: 1

Headphones

5 Pin Connector Headphones

Applicable Instruments: MEB-9100, MEB-9200/9300, MEB-9400

Part No: DR-531B-14

Qty: 1

Leg Wire Extensions

48-inch

Part No: C927

5/PKG

Ear Insert, Adult

Replacement Ear Inserts

Applicable Instruments: YE-102J

Part No: AAA-90080

Qty: 2/PKG

Ear Insert, Child

Replacement Ear Inserts

Application Instruments: YE-102J

Part No: AAA-90081

Qty: 2/PKG

Ear Insert, Set of 3 Sizes

Replacement Ear Inserts; includes 3 sizes in each package – SM, MED, LG

Applicable Instruments: YE-103J

Part No: NK-SVC-200

Qty: 1/PKG

Section

1

Electrode
Application
Supplies

2

EEG/PSG

3

EP/EMG/NCV