

<u>Manufacturer</u>	<u>Model</u>	<u>Interpretation</u>	<u>ECG Preview</u>	<u>12-lead Simultaneous Viewing</u>
Nihon Kohden	1250A	Yes	Yes	Yes
Nihon Kohden	1350A	Yes	Yes	Yes
Nihon Kohden	1550A	Yes	Yes	Yes
G. E.	MAC 1200	Yes	Yes	NO
Schiller	AT-1	Optional	NO	NO
Schiller	AT-2	Yes	NO	NO
Schiller	AT-2 LIGHT	Yes	NO	NO
Schiller	AT-2 PLUS	Yes	Yes	NO
Schiller	AT-101	Yes	Yes	NO
Schiller	AT-102	Yes	Yes	NO
Schiller	AT-10 PLUS	Yes	Yes	Yes
Welch Allyn	CP-100	Yes	Yes	NO
Welch Allyn	CP-200	Yes	Yes	NO
Burdick	Atria 3100	Yes	NO	NO
Burdick	Atria 6100	Yes	Yes	Yes
Cardioline	AR-600	Yes	NO	NO
Cardioline	AR-1200 view	Yes	Yes	NO
Cardioline	AR-1200adv	Yes	NO	NO
	Delta One			
Cardioline	PLUS	NO	Yes	NO
Cardioline	Delta Three	NO	Yes	NO
Cardioline	AR2100 view	Yes	NO	NO
Bionet	Model 2000	Yes	NO	NO
Bionet	Model 3000	Yes	Yes	NO
Fukuda Denshi	FCP7101	Optional	Yes	YES
Fukuda Denshi	FX7202	Optional	Yes	YES
Fukuda Denshi	FX7102	Optional	YES	YES
Fukuda Denshi	FX7402	Yes	Yes	YES

<u>Manufacturer</u>	<u>Model</u>	<u>Color Screen</u>	<u>8.5" X 11" Printout</u>	<u>ECG 3X4+1 Format</u>
Nihon Kohden	1250A	Yes	NO	Yes
Nihon Kohden	1350A	Yes	Yes	Yes
Nihon Kohden	1550A	Yes	Yes	Yes
G. E.	MAC 1200	NO	Yes	Yes
Schiller	AT-1	NO	NO	NO
Schiller	AT-2	NO	Yes	Yes
Schiller	AT-2 LIGHT	NO	Yes	Yes
Schiller	AT-2 PLUS	NO	Yes	Yes
Schiller	AT-101	NO	NO	NO
Schiller	AT-102	NO	Yes	Yes
Schiller	AT-10 PLUS	Yes	Yes	Yes
Welch Allyn	CP-100	NO	Yes	Yes
Welch Allyn	CP-200	Yes	Yes	Yes
Burdick	Atria 3100	NO	Yes	Yes
Burdick	Atria 6100	Yes	Yes	Yes
Cardioline	AR-600	No	NO	NO
Cardioline	AR-1200 view	NO	NO	NO
Cardioline	AR-1200adv	NO	NO	NO

Cardioline	Delta One PLUS	NO	NO	NO
Cardioline	Delta Three	NO	NO	NO
Cardioline	AR2100 view	No	Yes	Yes
Bionet	Model 2000	No	Yes	Yes
Bionet	Model 3000	NO	Yes	Yes
Fukuda Denshi	FCP7101	NO	NO	NO
Fukuda Denshi	FX7202	NO	NO	NO
Fukuda Denshi	FX7102	NO	NO	NO
Fukuda Denshi	FX7402	NO	Yes	YES

<u>Manufacturer</u>	<u>Model</u>	<u>Memory</u>	<u>Data Output</u>	<u>SD Card Software Upgradeable</u>
Nihon Kohden	1250A/1350A	3000 w/\$8 1-GB SD Card	Yes	Yes
Nihon Kohden	1550A	3000 W/ SD CARD	Yes	YES
G. E.	MAC 1200	Optional	Optional	NO
Schiller	AT-1	NO	NO	NO
Schiller	AT-2	NO	NO	NO
Schiller	AT-2 LIGHT	Optional	Yes	NO
Schiller	AT-2 PLUS	Optional	Yes	NO
Schiller	AT-101	Optional	Yes	NO
Schiller	AT-102	40	Yes	NO
Schiller	AT-10 PLUS	400	Yes	NO
Welch Allyn	CP-100	40	Yes	Yes
Welch Allyn	CP-200	40	Yes	Yes
Burdick	Atria 3100	40	Yes	NO
Burdick	Atria 6100	40	Yes	NO
Cardioline	AR-600	1	NO	NO
Cardioline	AR-1200 view	1	NO	NO
Cardioline	AR-1200adv	1	NO	NO
Cardioline	Delta One PLUS	1	NO	NO
Cardioline	Delta Three	1	NO	NO
Cardioline	AR2100 view	1	NO	NO
Bionet	Model 2000	1	Yes	NO
Bionet	Model 3000	1	Yes	NO
Fukuda Denshi	FCP7101	YES	Yes	PCMCI
Fukuda Denshi	FX7202	YES	Yes	PCMCI
Fukuda Denshi	FX7102	YES	YES	PCMCI
Fukuda Denshi	FX7402	Yes	Yes	PCMCI

<u>Manufacturer</u>	<u>Model</u>	<u>Warranty</u>	<u>Built-in ECG Waveform Simulator For User Training and Testing</u>	<u>Optional 15-lead</u>
Nihon Kohden	1250A	4 Year	Yes	NO
Nihon Kohden	1350A	4 Year	Yes	NO
Nihon Kohden	1550A	3 Year	YES	YES

G. E.	MAC 1200	3 Year	No	NO
Schiller	AT-1	3 Year	NO	NO
Schiller	AT-2	3 Year	NO	NO
Schiller	AT-2 LIGHT	3 Year	NO	NO
Schiller	AT-2 PLUS	3 Year	NO	NO
Schiller	AT-101	3 Year	NO	NO
Schiller	AT-102	3 Year	NO	NO
Schiller	AT-10 PLUS	3 Year	NO	NO
Welch Allyn	CP-100	3 Year	NO	NO
Welch Allyn	CP-200	3 Year	NO	NO
Burdick	Atria 3100	3 Year	NO	NO
Burdick	Atria 6100	3 Year	NO	NO
Cardioline	AR-600	2 Year	NO	NO
Cardioline	AR-1200 view	2 Year	NO	NO
Cardioline	AR-1200adv	2 Year	NO	NO
	Delta One			
Cardioline	PLUS	2 Year	NO	NO
Cardioline	Delta Three	2 Year	NO	NO
Cardioline	AR2100 view	2 Year	NO	NO
Bionet	Model 2000	3 Year	NO	NO
Bionet	Model 3000	3 Year	NO	NO
Fukuda Denshi	FCP7101	2 Years	NO	NO
Fukuda Denshi	FX7202	2 Years	NO	NO
Fukuda Denshi	FX7102	2 Years	NO	NO
Fukuda Denshi	FX7402	2 Years	NO	NO

Manufacturer	Model	Optional Stress Testing	12 inch Color Touch Screen	Internal Rechargeable Battery
Nihon Kohden	1250A	NO	NO	Yes
Nihon Kohden	1350A	NO	NO	Yes
Nihon Kohden	1550A	YES	YES	Yes
G. E.	MAC 1200	NO	NO	Yes
Schiller	AT-1	NO	NO	Yes
Schiller	AT-2	NO	NO	Yes
Schiller	AT-2 LIGHT	NO	NO	Yes
Schiller	AT-2 PLUS	NO	NO	Yes
Schiller	AT-101	NO	NO	Yes
Schiller	AT-102	YES	NO	Yes
Schiller	AT-10 PLUS	YES	NO	Yes
Welch Allyn	CP-100	NO	NO	Yes
Welch Allyn	CP-200	NO	NO	Yes
Burdick	Atria 3100	NO	NO	Yes
Burdick	Atria 6100	NO	NO	Yes
Cardioline	AR-600	NO	NO	Yes
Cardioline	AR-1200 view	NO	NO	Yes
Cardioline	AR-1200adv	NO	NO	Yes
	Delta One			
Cardioline	PLUS	NO	NO	Yes
Cardioline	Delta Three	NO	NO	Yes
Cardioline	AR2100 view	NO	NO	Yes
Bionet	Model 2000	NO	NO	Yes
Bionet	Model 3000	NO	NO	Yes
Fukuda Denshi	FCP7101	NO	NO	Optional

Fukuda Denshi	FX7202	NO	NO	Optional
Fukuda Denshi	FX7102	NO	NO	Optional
Fukuda Denshi	FX7402	Yes	NO	Optional

Manufacturer	Model	Wireless LAN Connectivity	Bar Coding	Analog Output
Nihon Kohden	1250A	YES	YES	Yes
Nihon Kohden	1350A	YES	YES	Yes
Nihon Kohden	1550A	YES	YES	Yes
G. E.	MAC 1200	Optional	Optional	NO
Schiller	AT-1	NO	NO	NO
Schiller	AT-2	NO	NO	NO
Schiller	AT-2 LIGHT	NO	NO	NO
Schiller	AT-2 PLUS	NO	NO	NO
Schiller	AT-101	NO	NO	NO
Schiller	AT-102	NO	NO	NO
Schiller	AT-10 PLUS	NO	NO	NO
Welch Allyn	CP-100	NO	NO	NO
Welch Allyn	CP-200	YES	NO	NO
Burdick	Atria 3100	NO	NO	NO
Burdick	Atria 6100	YES	Optional	NO
Cardioline	AR-600	NO	NO	NO
Cardioline	AR-1200 view	NO	NO	NO
Cardioline	AR-1200adv	NO	NO	NO
	Delta One			
Cardioline	PLUS	NO	NO	NO
Cardioline	Delta Three	NO	NO	NO
Cardioline	AR2100 view	NO	NO	NO
Bionet	Model 2000	NO	NO	NO
Bionet	Model 3000	NO	NO	NO
Fukuda Denshi	FCP7101	NO	NO	NO
Fukuda Denshi	FX7202	NO	NO	NO
Fukuda Denshi	FX7102	NO	NO	NO
Fukuda Denshi	FX7402	NO	NO	NO

Manufacturer	Model	Vectorcardiography	Apx. Street Price
Nihon Kohden	1250A	NO	\$1,982.00
Nihon Kohden	1350A	NO	\$2,722.00
Nihon Kohden	1550A	YES	\$ 5,995.00 (\$ 14,995.00 With STRESS incl. Treadmill)
G. E.	MAC 1200	NO	\$3,400.00
Schiller	AT-1	NO	\$1,250.00
Schiller	AT-2	NO	\$1,600.00
Schiller	AT-2 LIGHT	NO	\$2,200.00
Schiller	AT-2 PLUS	NO	\$2,900.00
Schiller	AT-101	NO	\$3,200.00
Schiller	AT-102	NO	\$5,400.00
Schiller	AT-10 PLUS	NO	\$7,000.00
Welch Allyn	CP-100	NO	\$4,900.00
Welch Allyn	CP-200	NO	\$5,400.00
Burdick	Atria 3100	NO	\$3,900.00
Burdick	Atria 6100	NO	\$4,900.00

Cardioline	AR-600	NO	\$1,199.00
Cardioline	AR-1200 view	NO	\$2,400.00
Cardioline	AR-1200adv	NO	\$2,200.00
Cardioline	Delta One PLUS	NO	\$1,399.00
Cardioline	Delta Three	NO	\$1,850.00
Cardioline	AR2100 view	NO	\$2,600.00
Bionet	Model 2000	NO	\$1,450.00
Bionet	Model 3000	NO	\$1,900.00
Fukuda Denshi	FCP7101	NO	\$2,495.00
Fukuda Denshi	FX7202	NO	\$3,000.00
Fukuda Denshi	FX7102	NO	\$1,995.00
Fukuda Denshi	FX7402	NO	\$3,850.00